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**Description**

curea P1 duo active is a sterile, hydro-active wound dressing. The surface is made of polypropylene fleece. The absorbent pad consists of cellulose composite with fluid-retaining sodium polyacrylate (SuperCore®) and a layer of activated carbon.

curea P1 duo active absorbs on both sides. P1 duo active absorbs exudate, blood and germs while maintaining a moist wound milieu even under compression.

The fluids are bound in the interior of the absorbent core. The absorption- and retention property of curea P1 duo active allows the wound a longer resting phase between the dressing changes thereby supporting the body's own wound healing. The product does not contain any pharmacological ingredients.

**Indication**

Chronic and low to heavy exuding wounds. Acute and iatrogenic wounds. Infected wounds or wounds vulnerable to infection. Malodorous wounds; Ulcerated malignant tumours.

**Contraindication**

Dry wounds (stagnation, risk of adhesion). Do not use for mucous membranes, eyes, sinews or bones, because these body parts might dry out and take permanent damage. Caution with cavities, because the wound pad expands when absorbing fluid.

**Application**

- curea P1 duo active offers two applications:
1. For odour control in wounds with clear exudate, put curea P1 duo active with the bright side against the wound;
  2. For decontamination put the dark side against the wound.

Make sure that curea P1 duo active is used under medical supervision and based on medical indication.

The physician treating the patient will define the intervals after which the dressing has to be changed based on the actual wound situation regarding exudation and infection. We recommend changing the dressing at least after 7 days and a frequent (best daily) check of

the wound area for signs of infection and of the dressing for depletion of capacity. Ensure that a sufficiently large area is available to attach the wound dressing on the skin.

**Adverse reactions**

No adverse reactions are known.

**Warnings**

- Do not cut or tear the wound pad. Polyacrylate might be spilled, contaminating the patient. In case polyacrylate is spilled, wash the patient with clean water and the wound with a wound irrigation solution.
- Dispose of damaged or expired products to avoid the risk of infection.
- The product is for single use only. Do not re-use to prevent the possible risk of infection.
- Do not re-sterilise.

**Other notices**

- The product has been sterilised with ethylene oxide.
- Latex-free

**Storage**

- Store in a dry place and protect from high humidity.
- Protect from direct sun light.
- Regard temperature limits.

**Disposal**

The wound pad can be disposed of with the household waste or in accordance with the local provisions. Special guidelines in medical institutions governing the disposal (e.g. for multi-resistant MRSA or VRE) must be complied with.

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