## MedCu Antimicrobial Wound Dressings with Copper-Oxide (AMWDs)

### Description

MedCu Antimicrobial Wound Dressings with Copper-Oxide (hereafter termed: AMWDs) are antimicrobial absorbent 2- or 3-layer sterile, soft, non-adherent, single use wound dressings consisting of the following:

- an internal absorbent cellulose and polyester layer containing copper oxide
- one or two external non-adherent nonwoven polypropylene layers containing copper oxide

The one or two external layers cover the internal layer from one or both sides, accordingly. The AMWD is available in several sizes, such as 10 cm x 10 cm, or 5 cm x 5 cm. Other sizes may be available. Consult individual device labels for exact dimensions. The wound dressing may or may not have an adhesive contour. The wound dressing is provided sterile using ETO sterilization method. Sterility is maintained as long as the pouch is not opened or breached. MedCu dressings absorb wound exudates, protect against contamination of the dressing, and serves as a barrier against external contamination reducing the risk of wound contamination.

### Indications

The MedCu Antimicrobial Wound Dressings with Copper-Oxide are indicated for the management of wounds and to provide an antimicrobial barrier. The dressings are applied topically and are in direct contact with the wound. The dressings are intended for indications as follows:

#### Over-The-Counter (OTC) use indications:

Local management of superficial wounds, minor burns, abrasions and lacerations (tears).

#### Prescription use indications (Rx-only):

Partial and full thickness wounds, pressure ulcers (stage I-IV), diabetic ulcers, venous stasis ulcers, arterial ulcers, 1<sup>st</sup> and 2<sup>nd</sup> degree burns, surgical wounds, vascular access or peripheral IV sites, orthopedic external pin sites, and acute wounds such as lacerations, abrasions, and skin tears.

### Contra-indications

Do not use on patients known to be hypersensitive to copper.

### For Rx-only use

#### Preparing the wound bed:

- Cleanse and dry peri-wound area thoroughly before dressing application as per facility protocol.
- If necessary, debride necrotic tissue before applying dressing.
- Hydrate wounds with low or no exudate as per facility protocol.
- If necessary, treat infection as per facility protocol.

### For both OTC and Rx-only

#### Applying the dressing:

- Choose the appropriate size dressing, peel-open the sterile pouch from the top side, and remove aseptically the wound dressing.
- If using the MedCu 3-layer wound dressings (orange on both sides), apply either side directly onto the wound. Cover with a secondary dressing and secure with fixation tape or a wrap.
- If using the MedCu 2-layer wound dressings (with one side orange) apply the orange side of the dressing to the wound site. Cover with a secondary dressing and secure with fixation tape or a wrap.
- If using the MedCu wound dressings with adhesive contour, no additional fixation, tape or wrap, are needed.
- The clinician should evaluate each wound and patient to determine best fixation method for the MedCu dressings.
- The clinician may choose any fixation method per institutional protocols.

### Frequency of dressing change:

- Frequency of wound dressing change may vary according to wound condition and clinician considerations.
- MedCu AMWDs may be left in place for up to 7 days.

### Dressing change and removal:


- MedCu AMWDs should be changed when signs of saturation of the dressing are noticeable or whenever clinician practice and facility protocol dictates.
- To change the non-adhesive MedCu dressings, remove the fixation bandages or tape and carefully lift the dressing off of the wound. Wet with sterilize saline, if needed, for easier removal of the wound dressing.
- To change the adhesive MedCu dressings, carefully peel the contour layer and gently lift the dressing off of the wound.

### Precautions

- Single use only • Do not resterilize • For external use only • Store at room temperature • Sterility is maintained as long as the sterile pouch is not breached
- Do not use if package is open or damaged • MedCu AMWDs may be used on infected wounds being managed in accordance with institutional clinical protocols for infection abatement as an adjunct to the standard treatment regimen • OTC use is allowed only for **superficial wounds, minor burns, abrasions and lacerations. Other wound types are Rx only.**

### Symbols interpretation

Symbols on the product package should be interpreted as follows:

	Expiration date (use by)		Single use – Do not re-use
	Consult operating instructions		Do not resterilize
	Do not use if sterile package compromised		Date of manufacture
	Sterilized using ETO		Manufacturer
	Store in dry conditions		Batch code / number
	Eu CE Symbol		Catalogue number

### If any additional information is required contact:

MedCu Technologies Ltd.  
Hasadnaot 10, Herzliya 4672837,  
Israel. [info@medcu.com](mailto:info@medcu.com)  
[www.medcu.com](http://www.medcu.com)

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