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## Patient Information Leaflet Revamil Wound Dressing

**The product**
Revamil Wound Dressing 5x5/8x8/10x20 cm honey-impregnated gauze for the treatment of wounds

**Product description**
Revamil Wound Dressing is a gauze impregnated with pure honey that has been produced under controlled conditions. The carefully controlled production process guarantees that Revamil Wound Dressing does not contain any traces of pesticides and that the same quality can be reproduced consistently. Revamil honey contains a high proportion of honey enzymes, such as glucose oxidase. The highly concentrated sugars cause some of the wound exudate to mix with Revamil gel. The resulting layer helps to strengthen the antibacterial protection provided by the dressing.

**Instructions for use**

- When required first clean the wound in the conventional way
- Open the packaging at the upper side by tearing open the seal
- Remove the impregnated dressing from the protective sheets and spread it out over the wound surface
- Attach the dressing in the normal way, e.g. by self-adhesive, porous tape
- Revamil Wound Dressing can also be covered by a secondary dressing that is attached in the conventional way
- After one to at most three days remove the dressing carefully from the wound and replace it with a new Revamil Wound Dressing

**Indications**
Revamil wound dressing can be used to treat different types of wounds such as bed sores and different kinds of ulcers, infected wounds, acute wounds, surgical oncological wounds, radiation induced oncological wounds and minor (first- and second-degree) burns.

**Contra indications**

- Do not use the product if the patient is known to be allergic or hypersensitive to honey. An allergic reaction or hypersensitivity can cause redness, irritation and possibly minor pain in the skin or wound.
- In the event of unexpected side effects after use, contact your doctor or the manufacturer. Side effects can cause inconvenience to the patient.

**Safety precautions**

- Do not use if honey is leaking from the inner packaging. Leakage can lead to contamination of the product with bacteria and to wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not use if the packaging (pouch) is not tightly sealed. Packaging that is not tightly sealed can lead to non-sterility of the wound dressing and to wound infection. Because of the high antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not apply this product for more than 30 days continuously. Treatment for more than 30 days continuously could lead to inconvenience to the patient, such as irritation and/or minor pain.
- Do not reuse an already used dressing. Reuse of a dressing can lead to loss of function of the honey, contamination of the product and/or contamination or infection of the wound.
- Do not sterilize an opened packaging. Resterilization of an opened packaging can lead to loss of function of the honey, contamination of the product and/or contamination or infection of the wound.
- Do not use the product after the expiration date as this can lead to lower performance of the product and to inconvenience for the patient.

**Storage conditions**
Store between 10 and 30°C

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